

Name: Ready-to-use ice pack	Code: FH1236	FH12 MEDICAL-HEALTH EQUIPMENT - FIRST AID
------------------------------------	---------------------	--

Description

Useful accessory to be used in the case of trauma or small bleeding; keeps the injured part cool for about 25/30 minutes depending on the outside temperature.

Disposable product.

Dimensions: 130 x 170 mm.

Medical equipment in class 2, complying to the requirements of Regulation CE 1907/2006 and subsequent modifications, included Regulation 552/09.



Category	Colours	Available Sizes In stock: U
Quantity/Carton 25	EC Standards	Important Information Suitable in case of trauma or minor bleeding. Read carefully the user sheet included. Make sure the device is in good condition before use. Once you become aware of a damage, you should immediately replace it. Keep in a dry place, far from heat sources and direct light.
Dexterity		Storage and cleaning