



Sir Safety System S.p.A. Single-member company Wholly-owned subsidiary of Sir Holding Srl Via dei Fornaciai, 9 - Santa Maria degli Angeli 06081 Assisi - Perugia - Italy Tel.+39 075.8043737 - Fax +39 075.8048686 VAT Nr. IT 03359340548 Viale Europa 74 20090 Cusago Milano Tel. +39 02.90394575 Fax +39 02.90394441

export@sirsafety.com www.sirsafety.com





Name: MEDIBAG 2 refill kit with sphygmoman Code: FH1225 FH12 MEDICAL-HEALTH EQUIPMENT - FIRST AID


Description

Supplies kit for Medibag 2 with sphygmomanometer.

Content

- no. 1 copy M.D. 388 15/7/03
- no. 5 pairs sterile gloves
- no. 1 mask with face shield
- no. 3 sterile saline solution, 500 ml
- no. 2 bottles of antiseptic solution, 500 ml, 10% Povidone-iodine
- no. 10 packs of sterile compressed gauze 10x10 cm
- no. 2 packs of sterile compressed gauze 18x40 cm
- no. 2 blankets 40x60 cm DIN 13152-BR
- no. 2 sterile forceps
- no. 1 tubular net gauze
- no. 1 pack of cotton wool
- no. 2 Plastosan boxes of 10 assorted patches
- no. 2 adhesive plaster rolls 5 m x 2.5 cm
- no. 1 Lister scissors 14.5 cm DIN 58279
- no. 3 tourniquets
- no. 2 ready-to-use Ice Packs
- no. 2 sanitary waste bags
- no. 1 digital clinical thermometer
- no. 1 sphygmomanometer with stethoscope
- no. 1 gauze bandage 3.5 m x 10 cm
- no. 1 First aid multilanguage instructions

First aid items required by It. M.D. 388 15/7/03 attachment 2 and It. L.D. 81 9/4/08 art. 45.



Category	Colours	Available Sizes
		Made to order: U In stock:
Quantity/Carton 1	EC Standards	Important Information
		Suitable for working environments, with three or more employees.
		Read carefully the user sheet included. Make sure the item isin good condition, undamaged before use. Once you become aware of a damage, you should immediately replace it.
		Keep in a dry place, far from heat sources.
Dexterity		Storage and cleaning